

# Ozone Therapy - Product Listing

## Ozone Therapy Device

**R9 400.00**



Ozone concentration: 99mg/L at 1/32 L/min oxygen; available for rectal and vaginal therapy and making of ozone drinking water; Built-in: cooling fan; ozone module; LED indicators with oxygen inlet connector and ozone outlet, connector OD 5mm; Driving voltage: 12 VDC (need 100-240 V AC)

## Mini Ozone Generator

**R10 350.00**



Destructor; Ozone concentration: 10 - 104mg/L; oxygen input: 1/32 - 1 lpm; available for rectal and vaginal therapy and making of ozone drinking water etc; Built-in: 3 levels ozone module; built-in ozone destructor, LED indicators, cooling fan, Connector: luer lock; Input voltage: 100-240 V AC 50/60 Hz to 12 V DC

## Double Quartz Ozone Generator

**R17 250.00**



Oxygen input: 1/32 - 1lpm; Built-in double quartz electrode ozone cell; built-in ozone destructor; Input voltage: 100-240 V AC 50/60 Hz to 12 V DC

## Non-return Valve

For air flow:  $\leq 5$ L/min

**R30.00**

High quality valve

**R90.00**



## Male Luer Connector

**R58.00**

5/32" fit for silicon hose with 4mm inner diameter



## Female Luer Connector

**R50.00**

5/32" fit for silicon hose with 4mm inner diameter



## Female Luer Connector Wall Mount

**R90.00**

Wall-mounted version with nut for locking





### Y-shape Plastic Connector

**R50.00**

Size: for 7x4mm silicon tubing



### 2-way Connector

**R172.00**

One end with 5/32" male luer,  
another end with 5/32 female luer



### 3-way Connector

**R195.00**

Two ends with 5/32" male luer,  
one end with 5/32 female luer



### Silicon Tube

Wall thickness: 1.5mm 7x4; Color: transparent; **R80.00**  
Material: silicon; ozone-resistant per metre

Wall thickness: 2mm 8x4; Color: transparent; **R95.00**  
Material silicon; ozone-resistant per metre



### Stone Diffuser

**R20.00**

Plastic mouth dia. 5.4mm; Material: sandstone



### White Stone Diffuser

**R40.00**

Plastic mouth dia. 5.4mm; Material: white sandstone



### Oxygen Tank Regulator

**R1 265.00**

(For use on MEDICAL oxygen tanks)

Regulate the flow of oxygen from the tank into the ozone generator; 1/32 l/m - 4 l/m; Maximum pressure - 3,000 PSI; Barb Outlet



### Oxygen Tank - 3 litre

**R2 800.00**

3 litre oxygen tank with medical tap  
Suitable for use with oxygen tank regulator



### Syringe

60 mL/CC; luer lock

**R20.00**

100 mL/CC - Glass; luer lock

**R680.00**



### Luer Lock Short Needles

**R15.00**

Luer lock connection; Needle length: 1/2 inch;  
Needle size: ID 1.2mm; OD 1.64mm



### Luer Lock Catheter

**R60.00**

Fr 12 size; with 2 eyes at the tip with luer connector  
single package for medical use



### Ozone Destructor for Ozone Therapy **R1 300.00**

Standard luer connector



### Ozone Bag for Limb Wound Treatment

**R900.00**

With bag adaptors (2 connectors & 2 bags)




### Ozone Bag for Insufflation

**R350.00**

With 3 chambers: 200ml, 400ml, 750ml with luer connector



**Insufflation kit**



**Olive oil breathing kit**



**Ozonide breathing kit**

**500 ml - R450.00**  
**1000 ml - R480.00**




**Ozonated water kit**

**Price on request**

*Prices are subject to exchange rate fluctuations*

October 2024



**Jochem Greiner - 074 1116658**  
**Susan Greiner - 083 500 8302**

**jochem@greinertrading.co.za**  
**www.greinertrading.co.za**

**Warranty and Disclaimer:**

All products have a one year warranty from date of purchase in which repairs or exchange of parts according to that which is deemed necessary by the company will be made free of charge. This warranty does not cover damage by the following:

- 1 Misuse, abuse or alteration to the product
- 2 Damage caused by accident or neglect
- 3 Servicing or repairs made by unauthorised persons

Neither Greiner Trading CC nor its Member, shall be held responsible for any incidental, special or consequential damages resulting from the use of any products. Treatment using ozone devices are not intended as a diagnosis, cure or as a substitute for professional medical advice.