# Ozone Therapy - Product Listing

# **Ozone Therapy Device**

R9 400.00



Ozone concentration: 99mg/L at 1/32 L/min oxygen; available for rectal and vaginal therapy and making of ozone drinking water; Built-in: cooling fan; ozone module; LED indicators with oxygen inlet connector and ozone outlet, connector OD 5mm; Driving voltage: 12 VDC (need 100-240 V AC)

#### Mini Ozone Generator

R10 350.00



Destructor; Ozone concentration: 10 - 104mg/L; oxygen input: 1/32 - 1 lpm; available for rectal and vaginal therapy and making of ozone drinking water etc; Built-in: 3 levels ozone module; built-in ozone destructor, LED indicators, cooling fan, Connector: luer lock; Input voltage: 100-240 V AC 50/60 Hz to 12 V DC

# **Double Quartz Ozone Generator**

R17 250.00



Oxygen input: 1/32 - 1lpm; Built-in double quartz electrode ozone cell; built-in ozone destructor; Input voltage: 100-240 V AC 50/60 Hz to 12 V DC

# -

#### Non-return Valve

For air flow: <=5L/min R30.00

High quality valve R90.00



#### Male Luer Connector

R58.00

5/32" fit for silicon hose with 4mm inner diameter



#### **Female Luer Connector**

R50.00

R90.00

5/32" fit for silicon hose with 4mm inner diameter



#### Female Luer Connector Wall Mount

Wall-mounted version with nut for locking



# **Y-shape Plastic Connector**

R50.00

Size: for 7x4mm silicon tubing



## 2-way Connector

R172.00

One end with 5/32" male luer, another end with 5/32 female luer



## **3-way Connector**

R195.00

Two ends with 5/32" male luer, one end with 5/32 female luer



#### Silicon Tube

Wall thickness: 1.5mm 7x4; Color: transparent;

Material: silicon; ozone-resistant

R95.00

R80.00

per metre

Wall thickness: 2mm 8x4; Color: transparent;

per metre

Material silicon; ozone-resistant



#### Stone Diffuser

R20.00

Plastic mouth dia. 5.4mm; Material: sandstone



#### White Stone Diffuser

R40.00

Plastic mouth dia. 5.4mm; Material: white sandstone



## **Oxygen Tank Regulator**

R1 265.00

(For use on MEDICAL oxygen tanks)

Regulate the flow of oxygen from the tank into the ozone generator; 1/32 l/m - 4 l/m; Maximum pressure - 3,000 PSI; Barb Outlet



### Oxygen Tank - 3 litre

R4 715.00

3 litre oxygen tank with medical tap Suitable for use with oxygen tank regulator



#### **Syringe**

60 mL/CC; luer lock **R20.00** 

100 mL/CC - Glass; luer lock **R680.00** 

Luer Lock Short Needles



Luer lock connection; Needle length: 1/2 inch;

Needle size: ID 1.2mm; OD 1.64mm



#### **Luer Lock Catheter**

R60.00

R15.00

Fr 12 size; with 2 eyes at the tip with luer connector single package for medical use



# Ozone Destructor for Ozone Therapy R1 300.00

Standard luer connector



# Ozone Bag for Limb Wound Treatment

R900.00

With bag adaptors (2 connectors & 2 bags)



#### Ozone Bag for Insufflation

R350.00

With 3 chambers: 200ml, 400ml, 750ml with luer connector



**Insufflation kit** 



Ozonide breathing kit

500 ml - R560.00 1000 ml - R600.00



Olive oil breathing kit




Ozonated water kit

**Price on request** 

#### Prices are subject to exchange rate fluctuations

November 2025





Jochem Greiner - 074 1116658 Susan Greiner - 083 500 8302

# jochem@greinertrading.co.za www.greinertrading.co.za

#### Warranty and Disclaimer:

All products have a one year warranty from date of purchase in which repairs or exchange of parts according to that which is deemed necessary by the company will be made free of charge. This warranty does not cover damage by the following:

- 1 Misuse, abuse or alteration to the product
- 2 Damage caused by accident or neglect
- 3 Servicing or repairs made by unauthorised persons

Neither Greiner Trading CC nor its Member, shall be held responsible for any incidental, special or consequential damages resulting from the use of any products. Treatment using ozone devices are not intended as a diagnosis, cure or as a substitute for professional medical advice.